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Burs are welded precisely to ensure best possible concentricity and stability. The blades are configured for specific applications to ensure that they are as gentle as possible when used in the patient's mouth.

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*The veneer should be cut with a diamond instrument first.*

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*The crown metal layer should be reduced gradually with the H34 carbide instrument.*

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

### When cutting crowns and bridges the following recommendations should be observed:

- 1 The veneer should be cut with a diamond instrument first.
- 2 Reduce the crown metal layer gradually, which means starting the procedure from top to bottom end several times. opt. speed 160,000 r.p.m.,  
Contact pressure: 0.5 N  
Please observe the working angle of 45 degrees!

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
*High cutting capacity blade design for separating bonded porcelain crowns made of low-fusing porcelain.*

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
### Features and advantages:

- Used for alloys faced with low-fusing porcelain
- The special blade design ensures high cutting capacity.
- The tip is designed for cutting and slices into the crown.
- One-piece production ensures high stability and reduces the risk of fracture.
- One instrument for two procedures - no need to change between diamond and carbide instruments.

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Caution: use spray coolant when preparing.

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
#### Crowns quickly separated through any material including zirconium oxide

- The 127  $\mu\text{m}$  grit size has a high cutting capacity that ensures all types of crowns can be quickly separated.
- The diamond grit cuts through gold, porcelain, non-precious metals and even zirconium oxide using one instrument.
- The disc facilitates separating bridges intraorally.
- With continuous cooling it swims over the surface of the gingiva, reducing the risk of damaging the gingival margin.

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**Important:** It is essential to cut a deep "V" into the crown. The cut is wider at the surface than deeper into the crown. Otherwise the disc can become stuck between the two sections of the crown, bend and/or break if the patient moves. Not suitable for the lifting of crowns (risk of fracture)!

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### The super sharp multifunctional cutting instrument.

#### For the areas of indication:

Crown cutting, removing fillings and preparing abutments.

The MF-Cut instruments have a special cutting geometry with a sharp-edged angle, which produces the extremely high cutting capacity on all types of alloy as well as amalgam and composites.

The end-cutting design of the instruments facilitates quick penetration of the material to be prepared.

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