

INgezim BTV CROM

R.12.BTV.K41

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

Ed.020217

INgezim BTV CROM R.12.BTV.K41

INgezim BTV CROM CROM is based on a direct Immunochromatography technique, which uses the VP7 protein of the Blue Tongue Virus to detect specific antibodies against such virus.

TECHNICAL BASIS

The device consists of a plastic casing with two windows:

- Window for sample addition: Contains the VP7 and a control protein, bound to coloured latex.
- Window for reading of results: Contains a test line (T) made of VP7 and a control line (C) made of an specific antibody to the control protein.

When the sample is added, if it contains specific antibodies to BTV, they will bind to the VP7 conjugated with latex and will flow along the membrane. The antigen-antibody-latex complex will bind to the antigen placed on the test zone (T) showing a pink line. The blue line in the control zone indicates the validation of the assay.

APPLICATION

Detection of antibodies specific to BTV VP7 in ruminant samples

DIAGNOSTIC SENSITIVITY

Study 1: 24 positive sera from animals vaccinated with different vaccines were analyzed and the results obtained were compared with INgezim BTV COMPAC. The correlation between both assays was 100%.

Study 2: 30 sera from animals (15 from cattle and 15 from lambs), experimentally infected and previously classified as positive by INgezim BTV COMPAC, were analyzed. The correlation with INgezim BTV COMPAC was 100%.

Study 3: A total of 389 positive field sera by INgezim BTV DR were analyzed. The results obtained indicated 97% of correlation between both assays. The 11 not matching sera showed OD values near the cut off by INGEZIM BTV DR.

DIAGNOSTIC SPECIFICITY

268 field sera from cattle from BTV free areas were analyzed. The specificity showed was 98%

ANALYTICAL SENSITIVITY

The following sera have been analyzed:

- OIE Reference Sera for serotypes 1, 2, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 20, 21, 22, 23, 24 (OIE Reference Laboratory, Pirbright).
- OIE Reference Serum serotype 19 (USDA).
- OIE Reference Sera serotype 3 (Spanish Reference Laboratory, Algete, Spain).
- 2 Animals infected and not vaccinated.
- 30 animals vaccinated (two doses: day 0 & 21).

The assay detects specific antibodies to the **24 serotypes of BTV studied**. On the other hand, the results obtained with sera from infected and vaccinated animals indicated that the assay detects antibodies **from day 5 post infection and from day 14 post vaccination**.

COMPOSITION OF THE KIT

- Chromatographic devices
- Vials with diluent

Product registered in Spain (1819-RD)

PRODUCT MANUFACTURED BY INGENASA

IT-73840 IT-73780 ISO 14001:2015 9191.INGE ISO 9001:2015 9175.ING2

SHELF LIFE: **24 months**
Stored at 4°C-25°C

Ed.020217