

INgezim IBR gE

R.12.IBE.K3

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

INgezim IBR gE

R.12.IBE.K3

INgezim IBR gE is an enzymatic assay based on a blocking ELISA technique, which uses a monoclonal antibody (MAb) specific of the Infectious Rhinotracheitis Virus (IBRV) gE protein, and inactivated antigens. This assay allows differentiating between infected animals and those vaccinated with a deleted vaccine.

TECHNICAL BASIS OF THE KIT

1. Plates are coated with inactivated IBR antigen. On these wells, serum samples are added and incubated.
2. If the samples contain specific antibodies to IBR gE protein, they will bind to the antigen.
3. When a MAb-HRPO specific to IBR gE protein is added, it will bind to the antigen in absence of blocking antibodies in the sample (negative animals or those vaccinated with deleted vaccines). If there are antibodies blocking the antigen (infected animals), the conjugate will not be able to bind to it. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

APPLICATION

Detection of specific antibodies to the Infectious Rhinotracheitis Virus in bovine serum and milk samples. The assay allows the differentiation between animals vaccinated with deleted vaccine and infected ones.

INTERPRETATION OF THE RESULTS

Two cut off are used for the results interpretation: positive and negative. Samples will be considered as **Positive** if their blocking % is equal to or higher than the positive cut off; **Negative** if their blocking % is lower than the negative cut off and **Doubtful** if the blocking % is between both values.

VALIDATION

1. Use of Reference Sera

- **OIE:** this assay has been evaluated with the OIE Reference sera EU1 (+), EU2 (+) and EU3 (-). Results obtained indicate that the assay fulfills the OIE requirements
- **FLI:** The assay has been evaluated with a panel of sera characterized by ELISA and SN (RC1, RC2, R3D, R31B, R32C) from the Friedrich Loeffler Institute. Correspondence with the expected results was 100%.


2. Use of milk

A total of 252 samples of milk tanks (80-200 samples) were analysed. These samples were from a herd vaccinated with deleted vaccine and with a remote history of not-deleted vaccine vaccination. All samples were assayed in parallel with INgezim IBR Compac (no differentiation between vaccinated and non-vaccinated animals).

All INgezim IBR gE positive samples were also positive by INgezim IBR Compac and there were not any sample negative by INgezim IBR Compac and positive by INgezim IBR gE.

3. Correlation with other commercial assays

- **Study 1 (internal):** 942 sera of animals from Spanish and other European countries farms were analyzed. Results were compared with the ones obtained by HerdCheck IBR gE Ab. INGEZIM IBR gE showed 93% sensitivity and 99% specificity in relation to HerdCheck IBR gE Ab and hence, 96% correlation.

- **Study 2 (external):** 41 Polish sera were analyzed by INgezim IBR gE and CVITEST Bovis IBRgE. Results obtained indicated 100% correlation between both assays.
- **Study 3 (external):** 80 Slovakian sera were analysed and results were compared with the ones obtained by HerdCheck IBR gE Ab. Both assays showed 95% correlation.

COMPOSITION OF THE KIT

- Microtitration plates of 96 well
- Vials with Positive Control
- Vials with Negative Control
- Vials with Peroxidase Conjugate
- Bottle with Washing Solution
- Bottle with diluent
- Bottle with stop solution
- Bottle with substrate (TMB) ready to use

REGISTRATION NUMBER 962 RD
PRODUCT MANUFACTURED BY INGENASA

SHELF LIFE: **18 months**
Stored at 2°C-8°C

Ed.020217