
**RIKTA®**

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# RIKTA® - 04/4

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EN

## ATTENTION!

- When purchasing the Device check the following:
  - the complete set of delivery (page 29);
  - compliance of serial numbers in this Operating Manual (page 43) with the numbers on the Device;
  - availability of sale date and stamp of the seller in the Warranty Card (page 43).
- **Do not switch on the Device if it was kept at temperatures below 0 °C!** The Device should be kept in conditions of room temperature for at least 3 hours.
- Before operating this Device, read this Manual carefully and keep it during the entire useful life of the Device.

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Dear customer!

We congratulate you on purchasing the RIKTA® – 04/4 device and we express confidence that it will help to maintain your health for many years!

**INVISIBLE LASER RADIATION!**  
**DO NOT STARE INTO THE BEAM AND DO NOT MAKE**  
**DIRECT OBSERVATIONS USING OPTICAL INSTRUMENTS!**  
**CLASS 1M LASER EQUIPMENT**

## TABLE OF CONTENTS

<b>RU</b> .....	3
<b>EN</b>	
TABLE OF CONTENTS .....	23
INTENDED USE .....	24
GENERAL INFORMATION .....	25
SAFETY PRECAUTIONS .....	26
DELIVERY SET .....	29
DEVICE .....	30
PREPARATION FOR OPERATION .....	32
DEVICE TURNING ON .....	32
DEVICE SERVICEABILITY CHECK (TESTING) .....	33
CONDUCTING THE THERAPY SESSION .....	34
SIMULTANEOUS USE OF TWO EMITTERS .....	35
DEVICE TURNING OFF .....	36
RADIATION POWER RECOMMENDATIONS .....	37
TECHNICAL SPECIFICATIONS .....	38
MAINTENANCE .....	39
POSSIBLE MALFUNCTIONS	
AND TROUBLESHOOTING .....	40
PACKAGE, STORAGE AND TRANSPORTATION .....	41
MANUFACTURER'S WARRANTY .....	42
ACCEPTANCE CERTIFICATE .....	43
GUARANTEE CARD .....	43

### DESIGNATIONS WHICH ARE USED IN OPERATING MANUAL

warning about the situation,  
in which it is possible  
to cause harm to health

mandatory action symbol

notes and tips

RIKTA®- 04/4 device is manufactured at the enterprise JSC «MILTA – PCP GIT», whose quality management system complies with the requirements of international standard ISO 13485:2016.

The Device is entered in the State register of medical devices of the Russian Federation; its registration certificate No. ФСР 2011/10170 has no expiration date.

## INTENDED USE

RIKTA®- 04/4 device is intended painless and noninvasive (without breaking the skin integument) treatment of wide range of human diseases. The Device can be used simultaneously with medicinal methods of treatment and other methods of therapy.

The Device is intended for use in polyclinics, clinics, hospitals, treatment rooms, as well as at home after consultation with a laser therapy specialist.

RIKTA®- 04/4 device provides simultaneous penetrative action on the organism of low-intensity laser radiation of the near-infrared range in combination with a constant magnetic field, visible red light and broadband infrared incoherent radiation (**photo-magnetic-infrared laser therapy**).

The treatment is carried out as **a contact method** (the emitter is setting on the skin), and **non-contact method** (when installing the emitter at a distance of 0.5–2 cm from the skin). In case of large size of the exposure area, treatment is performed by scanning-smooth movement of the emitter at a speed of 0.5–1.5 cm/sec.

## GENERAL INFORMATION

- The present operating manual (further – the Manual) is the document that certifies basic parameters and technical specifications of the magnetic-infrared laser therapy device RIKTA®- 04/4, guaranteed by the Manufacturer, and is intended for acquaintance with the operating rules of the Device. The manual is also a guide for technical service, transportation and storage.

**The Operating manual should be kept during all the Device service life.**

- It is necessary to read the Manual carefully before using the Device. Particular attention should be paid to compliance with the safety precautions listed on page 26, as well as to the list of medical contraindications to laser therapy.
- Methods of carrying out treatment procedures are stated in the “Methodical manual on application of the laser therapy device RIKTA®” (further, the Methodical manual), which is a part of the delivery set of each Device.
- The manufacturer reserves the right to make adjustments and amendments to this Manual and the technical documentation at any time without notice.
- The manufacturer reserves the right to make changes in the configuration and design of the Device, which will not affect its operation and basic technical specifications.

---

The Device and its components should be disposed in accordance with general applicable standards and local environment protection rules.

## SAFETY PRECAUTIONS

Using the Device is allowed only after acquaintance with this Manual and the Methodical manual.

### DON'T OPEN

- Never open the Device and never try to disassemble or modify its internal components. The Device has no components serviceable by the user. In case of malfunction, immediately stop using the Device and resort to the manufacturer or the service center.

### PLUG AND POWER CORD

- Do not connect the Device if the power cord or the emitter connecting cable is twisted or damaged.
- Do not insert the power cord plug into mains with wet hands to avoid electric shock.
- Insert the power cord plug in the mains socket completely until it is secure.
- Do not pull the cord when disconnecting the power cord, make sure to hold the plug when disconnecting the power cord.

### WHEN MALFUNCTIONS ARE OBSERVED

- In case of wear or damage of the cord or plug, as well as in case of strange odor and smoke, immediately turn off the Device, disconnect the power cord from mains and resort to the manufacturer or the service center.

### EXTERNAL EFFECTS

- It is prohibited to spray flammable substances in the vicinity of the Device.
- Avoid exposure of the Device to moisture and excessive heat.

### CHILDREN'S SAFETY

- Keep the Device away from children.

## WHEN USING THE DEVICE

- The Device falls within Class 1M laser equipment.


**When using the Device, it is recommended to avoid eyes exposure to direct or reflected radiation. Do not stare into the beam and do not make direct observations using optical instruments.**

- Medical personnel are recommended to wear goggles while operating the Device (available on request). It is allowed not to use goggles when carrying out a single treatment procedure on condition that the Device aperture is tightly pressed to a patient's body.  
If the treatment is carried out in a patient's face area with the help of light guide nozzles, the goggles are also recommended for the patient.
- Start the laser therapy procedure only after putting the Device aperture on the treatment zone.
- At the time of treatment, strictly follow the guidelines for setting the emitter power (see page 37).

## DEVICE CARE

- Public Health Regulations on Arrangement and Operation of Lasers No. 5804-91 should be observed in the case the Device is placed and operated in medical facilities.
- Before and after a treatment procedure, make sure to wipe the Device with a cloth of gauze or coarse calico soaked in a treatment tool disinfectant based on 3% hydrogen peroxide solution. Avoid Device contact with alcohol-containing media!
- The Device disinfection may be carried out only when the Device is turned off and disconnected from mains.

## OPERATION CONDITIONS

The Device is designed for operation at a temperature from +10 to +35 °C, atmospheric pressure from 630 to 800 mm Hg, and relative humidity not more than 80%.

## MEDICAL CONTRAINDICATIONS TO LASER THERAPY

The contraindications to the Device home use include diseases accompanied with apparent structural irreversible changes. The contraindications can also reasonably include situations where an emergency medical care is required. If there are associated diseases (endocrine, autonomic disturbances) it is recommended to consult the doctor for the advisability of laser therapy treatment.

It is **STRICTLY FORBIDDEN** to carry out laser therapy at home in presence of the following medical contraindications\*:

- individual incompatibility of influence factors of the Device;
- general severe condition of the patient;
- blood diseases with the syndrome of bleeding;
- thrombocytopenia below 60,000 per  $\mu\text{l}$ ;
- acute feverish conditions;
- active tuberculosis;
- decompensated condition of cardiovascular, excretory, respiratory and endocrine systems;
- deep vein thrombosis;
- malignant neoplasm;
- benign neoplasms with a tendency to grow and become malignant;
- presence of implanted pacemaker;
- all terms of pregnancy;
- period 2 days before and during menstruation.

In case of exposition on the heart area (front and side surfaces of the thoracic), the frequency should be 5 Hz strict.

\* See Methodical manual.

## DELIVERY SET

No	Description	Quantity
<b>1</b>	Electronic Control Unit	1 pc
<b>2</b>	Emitter T1-04	1 pc
<b>3</b>	Operating manual	1 pc
<b>4</b>	Methodical manual for using RIKTA®- 04/4 device	1 pc
<b>5</b>	Device storage bag	1 pc

Optional accessories at request:

- **laser-safety goggles**

### Optional emitters:

- **emitter T1-04** (*additional*)  
The second T1 emitter reduces the procedure time and increases efficiency due to simultaneous effect of two radiation sources.
- **emitter T2-04**  
Similar to T1 but features higher laser radiation power ( $20^{+2}_{-4}$  W).
- **emitter T1E-04**  
Similar to T1 but features an additional treatment factor – transcutaneous electrical nerve stimulation (TENS).
- **emitter T2E-04**  
Similar to T1E but features higher laser radiation power ( $20^{+2}_{-4}$  W).
- **emitter Douche 1-04**  
Larger exposure area and higher power (30 W) for treatment of large surfaces.
- **emitter Douche 1 (50)-04**  
Similar to Douche 1 but features higher laser radiation power (50 W).
- **emitter Douche 2-04**  
Indicated for treatment of diseases in hairy part of the head.


### Stes of optical light-guide nozzles:

- **KON-1**  
used for endocavitary exposure in gynecology, otorhinolaryngology, stomatology, urology, reflexotherapy, cosmetology, etc.
- **KON-3**  
used in stomatology.
- **KON-Г**  
used in gynecology.

## DEVICE

### RIKTA®-04/4

Basic configuration:  
Electronic control unit  
and emitter T1-04

**Figure legend:**

- |  |   |
|--|---|
| <ul style="list-style-type: none"> <li>1 – <b>Electronic control unit</b></li> <li>2 – manufacturer's trademark</li> <li>3 – laser radiation warning sign</li> <li>4 – important User Manual information</li> <li>5 – power ON/OFF switch</li> <li>6 – power supply cord</li> <li>7 – 1 and 2 emitter connectors</li> <li>8 – LED indicators:<br/>LASER – laser radiation;<br/>IRD – infrared radiation;<br/>RED – red radiation</li> <li>9 – FREQUENCY – button to set pulse frequency and set frequency LEDs</li> <li>10 – TIME – button to set procedure time and procedure time LEDs</li> <li>11 – POWER LASER – button to set laser radiation power and set laser radiation power LEDs</li> <li>12 – POWER IRD – button to set IR radiation power and set IR radiation power LEDs</li> <li>13 – START/STOP – treatment procedure START/STOP button</li> <li>14 – SESSION – emitter operating status indicators</li> </ul> | <ul style="list-style-type: none"> <li>15 – <b>Emitter T1</b></li> <li>16 – emitting port with thread to attach light-guide nozzles</li> <li>17 – nameplate</li> <li>18 – emitter power supply unit connector</li> <li>19 – treatment procedure START/STOP button</li> <li>20 – laser radiation warning sign</li> </ul> |
|--|---|

**Nameplate at the lower space of the housing of the electronic control unit**

- L1 – manufacturer's trademark
- L2 – laser radiation warning sign
- L3 – laser safety class
- L4 – electrical safety, protection class II
- L5 – electrical safety, type BF
- L6 – manufacturer
- L7 – specification number
- L8 – supply network parameters

## PREPARATION FOR OPERATION

After staying at low temperature the Device shall be kept at a room temperature in the range of  $+20 \pm 5^{\circ}\text{C}$  for at least 3 hours before switching it on.

### DEVICE TURNING ON

- Before connecting the Device make sure that the power ON/OFF switch [5] is in the OFF or “0” position and the emitter is connected to “1” or “2” [7] connector of the electronic control unit [1].
- Insert the power cord plug [6] in the mains socket and set the power ON/OFF switch [5] to the ON or “1” position.

At this point:

the following LEDs light up

FREQUENCY [9] – 5 Hz, TIME [10] – 5 min,

POWER LASER [11] – 100%, POWER IRD [12] – 100%;

an audible alarm sounds.

Prior to operation, it is recommended to test the Device (see «DEVICE SERVICEABILITY CHECK (TESTING)» on page 33).

## DEVICE SERVICEABILITY CHECK (TESTING)

The Device features the test mode to check emissions at the emitter aperture. To start the testing procedure:

- Apply the emitter [15] aperture [16] to any flat bright surface (for example, to a clean white paper sheet or palm), press the START/STOP button [19] on the emitter and hold it for 3–4 sec.

The Device switches to the self-test mode: a short audible alarm sounds, LEDs LASER, IRD, RED [8] light up one by one. The testing procedure lasts several seconds.

If, upon the test completion, a short audible alarm sounds again and LEDs LASER, IRD, RED [8] light up all together, then **the test is successful**. The Device is ready for operation.

If, upon the test completion, any of the LEDs [8] blinks while an intermittent audible alarm sounds, the emitter is **FAULTY** and cannot be further used.

Turn off the power (set the supply voltage switch [5] to the OFF or “0” position), disconnect the power cord plug from the mains socket.

Consult the manufacturer or service center.

## CONDUCTING THE THERAPY SESSION

- Based on the doctor's prescription or in accordance with the guidelines set out in the Guidance Manual, given the recommendations on page 37 hereof, use:
  - button FREQUENCY [9] to set the desired pulse frequency,
  - button TIME [10] to set the desired procedure time,
  - buttons POWER LASER [11] and POWER IRD [12] to set the desired laser and infrared radiation power.LEDs FREQUENCY, TIME, POWER LASER, POWER IRD will light up according to the selected settings.
- Apply the emitter to the zone to be treated.
- Shorty push the START/STOP button on the electronic control unit [13] or on the emitter [19].

An audible alarm sounds, LED SESSION [14] starts blinking:

*green* if the emitter is connected to the connector "1",

*red* if the emitter is connected to the connector "2",

*alternatively red and green* if both emitters are connected.

During the treatment procedure, LEDs TIME [10] will successively light up to display the remaining procedure time.

Upon the procedure completion an audible alarm sounds, LED SESSION [14] is off. LEDs TIME [10] will display the set procedure time again.

**Note:**

*The treatment procedure can be stopped at any time by pressing the START/STOP button on the electronic control unit [13] or on the emitter [19].*

## SIMULTANEOUS USE OF TWO EMITTERS

The Device supports simultaneous use of two emitters. The convenience of two emitters consists in that different exposure zones can be treated simultaneously with the same exposure parameters to reduce the total procedure time.

**Note:**

*For the list of emitter models compatible with the Device refer to page 29, “OPTIONAL EMITTERS”.*

Before connecting the Device to the mains, connect the emitters to “1” and “2” connectors. Then, proceed according to the description in section “PREPARATION FOR OPERATION” (page 32).

**Note:**

*When testing the Device (section “DEVICE SERVICEABILITY CHECK (TESTING)”, page 33), test each emitter separately.*

Operation with two emitters is similar to the operation with one emitter, as described in section “CONDUCTING THE THERAPY SESSION” (page 34) but for some peculiarities:

- START/STOP button [13] on the Device electronic control unit starts/stops both emitters at once.
- START/STOP button [19] on the emitter starts/stops only that emitter.
- LED SESSION [14] on the Device electronic control unit with concurrent operation of two emitters will be blinking alternatively red and green.

## DEVICE TURNING OFF

To turn off the Device set power ON/OFF switch [5] to the OFF or “0” position. All LEDs on the electronic control unit go out.

Upon completion of operation, disconnect the power cord plug [6] from the mains socket.

## RADIATION POWER RECOMMENDATIONS

 For children under 14, fragile patients after heart attack or cerebral accident, people above 70 and disabled persons, **use only low power emitter T1-04**, as listed in the table.

Recommended power of laser and infrared radiation depends on the patient's age and condition.

Age	T1-04 emitter laser and infrared radiation power
Newborn babies and children under 1 year	25% (Distant technique only!)
1-8 years	25%
8-14 years	50%
<b>Children above 14 and adults</b>	<b>100%</b>
Fragile patients after heart attack or cerebral accident	50%
People above 70	50%
People above 80 and disabled persons	25%

## TECHNICAL SPECIFICATIONS\*

### Emitter T1:

Radiation wavelength,  $\mu\text{m}$ :

- laser infrared radiation ..... 0,80...0,91
- broadband infrared radiation..... 0,86...0,96
- visible red radiation ..... 0,60...0,67

Laser radiation pulse power, W .....  $16^{+2}_{-4}$

Broadband infrared radiation average power, mW .....  $60 \pm 30$

Red radiation average power, mW .....  $7^{+3}_{-5}$

Magnetic field induction, mT .....  $35 \pm 10$

Device aperture area,  $\text{cm}^2$  .....  $4 \pm 0,4$

Laser and broadband infrared radiation pulse repetition frequency, Hz:

- fixed .....  $(5, 50, 1000) \pm 5\%$
- variable in the range..... 250...1

Red light pulse repetition frequency, Hz..... 2

Exposure time, min. ....  $(1, 2, 5, 10) \pm 5\%$

The Device power supply, AC network:

- voltage, V ..... 50
- frequency, Hz ..... 20

Electrical safety ..... protection class II

Applied part ..... type BF

Laser safety by IEC 60825-1 .....class 1M

The Device complies with IEC 60601-1-2:2014 standard on electromagnetic compatibility.

Control unit dimensions, mm, max ..... 245x220x95

Net weight, kg, no more ..... 1,5

Service life, years ..... 5

\* For reference.

## MAINTENANCE

- To ensure reliable operation of the Device, perform its maintenance in a timely manner, guided by the table 1.
- Maintenance shall be carried out before each use of the Device and shall include the following steps:

Table 1

Job content	Requirements
A. Visual inspection without using special tools	Absence of damages to the control unit, emitter, connecting cable and power cord
B. Cleaning and disinfection of the emitter in accordance with the requirements set out in section “SAFETY PRECAUTIONS”, page 26	
C. Check the following:	
1. Proper and tight connection of the power cord	The power cord should be attached to the plug and to the Device housing in such a way as to prevent the cable from twisting
2. Lit LEDs after the supply voltage switch is set to the ON position	The following LEDs should be lit: FREQUENCY – 5 Hz, TIME – 5 min, POWER LASER – 100%, POWER IR – 100%
3. Radiation at the emitter outlet (test the Device, as described on page 33)	Upon the test completion a short audible alarm sounds. LEDs LASER, IRD, RED [8] light up and go out all together

- In case of detection of noncompliance of the Device with the requirements specified in table 1, further operation of Device is not permitted, the Device is to be repaired.

## POSSIBLE MALFUNCTIONS AND TROUBLESHOOTING

The list of possible malfunctions, causes and troubleshooting are listed in table 2.

*Table 2*

Malfunctions	Causes	Troubleshooting
No LEDs light up after the supply voltage switch is on	1. No voltage in the mains socket	Check if there is voltage supplied to the socket
	2. Power cord fault	Repairs are carried out by the Manufacturer
	3. Device fault	Repairs are carried out by the Manufacturer
A LED does not light up at successive pushing of buttons FREQUENCY [9], TIME [10], POWER LASER [11], POWER IRD [12]	Trouble of normal operation of the Device	Repairs are carried out by the Manufacturer
On starting the treatment procedure ("CONDUCTING THE THERAPY SESSION", page 34) LED SESSION [14] is not blinking	Trouble of normal operation of the Device	Repairs are carried out by the Manufacturer

## PACKAGE, STORAGE AND TRANSPORTATION

- Each Device with accessories is supplied in the manufacturer's package.
- The Device should be stored in the Manufacturer's package in the heated room at temperature in the range of +10...+35 °C and relative humidity of no more than 80% (at +25 °C). The storage room should be free of dust, mercury vapor, acids and alkalis provoking corrosion.
- The Device should be transported by all types of the closed transport in Manufacturer's package. It is necessary to ensure a stable position of the package during transportation.

## MANUFACTURER'S WARRANTY

- The Manufacturer guarantees that the Device parameters correspond to the requirements of the Technical Specifications TU 9444-007-18537172-01, within 24 months from the date of sale of the Device to the consumer (if the consumer complies with the conditions of use, transportation and storage). Within the aforesaid period, the manufacturer provides free repairs of the Device.
- The Manufacturer is not responsible for and does not make free repairs of malfunctions caused by a fault of the consumer (for example, due to careless handling, transportation, non-compliance with the rules of use, as well as in case of breaking of the seal).
- This warranty does not apply to any of the following:
  - absence of the present Manual;
  - absence of the Warranty card;
  - mechanical damage including the damage during transportation;
  - the Device failure because of foreign objects getting inside it (nails, clips, etc.) and any liquids.
- The warranty also becomes void if during the warranty period a repair of the faulty Device was carried out by unauthorized persons and organizations.
- Upon the warranty period expiration, the consumer shall repair the Device at own expense.
- In the absence of the sales date and the seller's stamp in the Warranty card, the warranty period shall be calculated from the date of Manufacturing of the Device.
- All issues related to repairs and maintenance should be directed to the Manufacturer or to the technical centers listed on the website [www.rikta.ru](http://www.rikta.ru) or in a separate warranty booklet (if available).

After the end of service life, the Manufacturer is not responsible for the safety and compliance of the product with the functional characteristics stated in this Manual.

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