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The **Foot Aging Syndrome** is a set of conditions that, on the one hand, results in structural and functional foot changes connected with age and, on the other, on the localization of systemic chronic conditions that often (but not only) affect the elderly and which may influence mobility and quality of life. Among these conditions, we cite pain, minor or severe deformities, traumas or chronic diseases such as the **Diabetic Foot**. All this leads to restricted movement and loss of self-esteem, which, in turn, negatively affects social life. **Optima Molliter** carefully studies methods, materials and technologies and recreates the condition for ideal movement, providing a menu of technologies and methods that allow patients to make their choice. Today, our **Wound Care program** and its **Algorithm** help Doctors, Podiatrists and Orthopedic Technicians across the world to provide effective treatments, thus reducing costs and healing time.

**Our goal is to make people free
to move in their world!**

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“The **DFA** (diabetic foot offloading algorithm) is a “lingua franca” that can work in Shanghai and Abu Dhabi, and all the way to Los Angeles, Houston or Rome”. It is inspired by the fundamental principles and scientific elements of international guidelines. The **flexibility of this algorithm** is that it can mix and combine technologies and be easily understood by everybody, thus providing appropriate therapeutic guide to doctors, podiatrists or healthcare operators who take care of patient presenting the **risk, high risk and very high risk of diabetic foot**”.

*Prof. David G. Armstrong,
DMP, MD, PhD,
Professor of Surgery and Director,
Southwestern Academy Limb Salvage Alliance (SALSA)
Keck School of Medicine at University of Southern California*

Diabetic Foot Off-loading Algorithm[®]

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
Risk 0 Foot at no risk

No LOPS¹, no PAD² no Deformity


Risk 1 Foot at risk / Risk 2 Foot at high risk

LOPS¹ and/or PAD² and/or foot deformity and/or pre ulcerative signs

Risk 3 Foot at very high risk


LOPS¹ and/or PAD² + Severe deformity and/or history of ulceration and/or minor/major amputation and/or end-stage renal disease


ULCERATION: Plantar forefoot or midfoot

Neuropathic plantar forefoot or midfoot Ulcer

Recommendation Ia - Ib*

Recommendation II*

Recom. III*

Recommendation IVa - IVb*

Don't use standard or therapeutic shoes (IVa)

OPTIMA Free + Monolayer insole or Puzzle Insole[®] + felted foam

Relief Dual[®] + PegAssist[®] Insole or PegContour[®] Insole + felted foam

OPTIMA Europa + Monolayer insole or Puzzle Insole[®] + felted foam

WCS[®] + WCS[®] Insole + felted foam


AllRound Shoe[®] + PegContour[®] Insole + felted foam


Recommendation V* and Recommendation VI*

In case of neuropathic metatarsal head or digital ulcer consider prophylactic surgery and curative surgery + post operative therapeutic walkers or shoes

Ischemic and/or Infected plantar Ulcer Mild/Moderate/Severe

Recommendation VII a* Mild Infection or Mild Ischemia

Recommendation VII b* Mild Infection and Mild Ischemia or Moderate Infection or Moderate Ischemia

Recommendation VII c* Moderate Infection and Moderate Ischemia or Severe Infection or Severe Ischemia

Only after Ischemia and Infection addressing


ULCERATION: Heel and NON Plantar / MINOR AMPUTATION / CHARCOT FOOT

Neuropathic plantar Heel Ulceration


Recommendation VIII*


Neuropathic NON plantar foot Ulcer

Minor amputation management

Charcot foot

Essential Bibliography: A. Piaggini et al. Foot & Ankle International, April 15, 2016, vol. 37, Issue 8, 2016 | A. Piaggini, Diabetes Care, March 2007, vol. 30, no. 3 586-590 | J. Donnelly, JWC, vol. 20, Number 7, July 2011 | R. Dahmen, Diabetes Care 24: 705-709, 2001 | E. Faglia, G. Clerici, Foot Ankle Int., 2013 Feb; 34 (2): 222-7 | D. G. Armstrong, Diabetes Care 24: 1019-1022, 2001 | C. Caravaggi, E. Faglia, Diabetes Care 23: 1746-1751, 2000 | D. Armstrong, Diabetes Care, March 2005, vol. 28, no. 3 551-554 | Stephanie C. Wu, Diabetes Care 31: 2118-2119, 2008 | Ira A. Katz, Diabetes Care, March 2005, vol. 28, no. 3 555-559

¹ LOPS: Loss of Protective Sensation ² PAD: Peripheral artery disease

* Internat. Guidelines of Diabetic Foot: <https://iwgdfguidelines.org/guidelines/guidelines/>

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THE ONLY DEVICE WHICH HAS PROVED ITS EFFICACY IN THE HEALING OF THE DIABETIC PLANTAR LESIONS GRADE IA and IIA

(according to the Texas University classification)

Emerging Treatments and Technologies

ORIGINAL ARTICLE

An Off-the-Shelf Instant Contact Casting Device for the Management of Diabetic Foot Ulcers

A randomized prospective trial versus traditional fiberglass cast

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OBJECTIVE— This study was designed to test the safety, effectiveness, and costs of off-loading with a novel, off-the-shelf irremovable device in the management of diabetic foot ulceration (DFU).

RESEARCH DESIGN AND METHODS— We prospectively evaluated off-loading of neuropathic plantar ulcers in 40 diabetic outpatients attending our diabetic foot clinic and compared healing rates at the 12-week follow-up, number and severity of adverse events, healing time, costs and applicability of the device, and patients' satisfaction between those randomly assigned to total contact casting (TCC; group A) or to the Optima Diab walker (group B). Deep or infected ulcers were excluded.

RESULTS— No difference between groups A and B was observed in healing rates at 12 weeks (95 vs. 85%), healing time (6.5 4.4 vs. 6.7 3.4 weeks), and number of adverse events (six versus four). Treatment was significantly less expensive in group B, which showed a mean reduction of costs of 78% compared with group A (P 0.001). Practicability was more favorable in group B, with a reduction of 77 and 58% of the time required for application and removal of the devices, respectively (P 0.001). Patients' satisfaction with the treatment was higher in group B (P 0.01).

CONCLUSIONS— The Optima Diab walker is as safe and effective as TCC in the management of DFU, but its lower costs and better applicability may be of help in spreading the practice of off-loading among the centers that manage the diabetic foot.

Diabetes Care 30:586–590, 2007

Comparison of Removable and Irremovable Walking Boot to Total Contact Casting in Off-loading the Neuropathic Diabetic Foot Ulceration

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ABSTRACT

Background: Despite its efficacy in healing neuropathic diabetic foot ulcers (DFUs), total contact cast (TCC) is often underused because of technical limitations and poor patient acceptance. We compared TCC to irremovable and removable commercially available walking boots for DFU offloading.

Methods: We prospectively studied 60 patients with DFUs, randomly assigned to 3 different offloading modalities: TCC (group A), walking boot rendered irremovable (i-RWD; group B), and removable walking boot (RWD; group C). Patients were followed up weekly for 90 days or up to complete re-epithelization; ulcer survival, healing time, and ulcer size reduction (USR) were considered for efficacy, whereas number of adverse events was considered for safety. Patients' acceptance and costs were also evaluated.

Results: Mean healing time in the 3 groups did not differ (P = .5579), and survival analysis showed no difference between the groups (logrank test P = .8270). USR from baseline to the end of follow-up was significant (P < .01) in all groups without differences between the groups. Seven patients in group A (35%), 2 in group B (10%), and 1 in group C (5%) (Fisher exact test P = .0436 group A vs group C) reported nonsevere adverse events. Patients' acceptance and costs were significantly better in group C (P < .05).

Conclusions: Our results suggest that a walking boot was as effective and safe as TCC in offloading the neuropathic DFUs, irrespective of removability. The better acceptability and lesser costs of a removable device may actually extend the possibilities of providing adequate off-loading.

Level of Evidence: Level II, prospective comparative study.

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