
Medical Imaging Workstation

Operation Manual

Thanks for purchasing ezcap292 HD Medical Imaging Workstation. The operation manual guides you how to set up and operate.

Features/ports and buttons


Front panel:

1. HOST – Plug a USB flash (or USB HDD) to save recorded video files or update firmware.
2. IR – Receive the command signal from remote control.
3. VIDEO – Press to start recording video, press again to stop recording.
4. PHOTO – Press to snapshot a photo. **(No effect during recording video)**
5. MENU – Press to enter the setup menu.
6. EXIT – Quit current function.
7. PLAYBACK – Press to switch from Preview record mode to Playback mode.
8. PLAY/PAUSE – Press to play video file in Playback mode, and press again to pause.
9. DELETE – Delete file in Playback mode.
10. HDD/USB Switch – Press to switch between HDD and USB HOST.
11. POWER – Power on when connect power supply, press to power off, and press again to power on again. **(Device will remember all the settings before power off when you power on again. No effect during recording.)**

12. OK/UP/DOWN/LEFT/RIGHT – For menu operation. UP/DOWN also for volume up and down, and LEFT/RIGHT also for microphone volume up and down.

Back plate:

13. FOOT SWITCH – Foot pedal ports, PHOTO for snapshot, and VIDEO for recording video.

14. LINE OUT/LINE IN – Audio line output / Audio line input.

15. IR – Connect IR receiver cable and fix the receiver at somewhere for easy use of remote control.

16. MIC – Microphone input.

17. SDI OUT / IN – SDI output / input.

18. YPbPr IN – Component video (YPbPr) input.

19. CVBS OUT / IN – CVBS output / input.

20. HDMI IN – HDMI input.

21. DVI IN – DVI input.

22. VGA IN – VGA input.

23. PC – Connect to USB 3.0 port of computer with USB cable, to capture video with medical software.

24. HDD – Connect 2.5 inch HDD.

25. HDMI OUT – HDMI output.

26. RS232 – Extension port.

27. READ HDD – Connect to computer to read the data in 2.5 inch HDD with mini USB cable.

* **Only in USB HOST mode. And cannot read HDD in HDD mode.**

28. POWER – Connect DC 12V3A power supply.

Remote Control:

1. POWER – Power on and off.

2. PHOTO - Press to snapshot a photo. **(No effect during recording video)**

3. HDMI/YPbPr/AV/SDI/DVI/VGA – Select video source.

4. USB – Select USB HOST to be storage.

5. HDD – Select HDD to be storage.

6. REC/STOP – Record video / Stop recording.

7. OK/UP/DOWN/LEFT/RIGHT –For menu operation.
UP/DOWN also for volume up and down, and LEFT/RIGHT also for microphone volume up and down.

8. MENU - Press to enter the setup menu.

9. EXIT - Quit current function.


10. WATERMARK – Press to add time watermark with recorded file on top left, and press again to cancel watermark.

11. PLAYBACK - Press to switch from preview record mode to playback mode.

12. SCH REC – Set up schedule recording.

13. PLAY/PAUSE/LAST/NEXT/FAST BACKWARD/FAST FORWARD/DELETE – Function keys in Playback mode.
DELETE button can delete the selected file.

14. 720P/1080P/480P – Set up the resolution of video recording and HDMI bypass.

Connections

For connection, package only includes Power supply, Remote control, IR receiver cable, USB 3.0 cable, Mini USB cable. Please prepare yourself for other cables and accessories.

1, Connect DC12V power supply to POWER port of back plate, the green LED next to POWER port turns on, and the green LED next to POWER button of front panel turns on too.

2, Connect video source, and you have following selections: CVBS, YPbPr, SDI, HDMI, DVI or

VGA.

Please connect one of them only at the same time, and don't connect them together. The video source that has been connected first is effect if you connect two kinds of video sources or above at the same time. The device has the memory function, it will choose the video source that has been selected before device powered off. If you do not connect this video source, device will search other connected video sources. You also can select the video source by remote control.

3, Connect audio source, and you have following selections: audio line input, Microphone input.

4, Connect video output to display, and you have following selections: CVBS output, SDI output, HDMI output.

It's not necessary because device includes screen.

HDMI output is effect for all the video sources input, but CVBS output is effect for CVBS input only, and SDI output is effect for SDI input only.

5, Connect LINE OUT to speaker or headset.

Please connect this if you need to hear sound. And don't need connect it if video source is HDMI input or SDI input, and they include sound, and connect display with speaker.

6, Connect storage device, and you can choose USB flash drive or USB HDD to connect USB HOST port, or choose 2.5 inch HDD to connect HDD port. Please press USB/HDD switch to choose the storage mode, and the LED next to the mode will turn on.

Device supports NTFS, FAT and FAT32 format, and doesn't support exFAT format.

7, Connect foot pedal to FOOT SWITCH port. PHOTO port for snapshot photo, and VIDEO port for start/stop recording video.

8, Connect IR receiver cable to IR port of back plate, and fix IR receiver at somewhere. It's not necessary if you don't need.

9, Connect USB 3.0 cable to PC port of back plate and free USB 3.0 port of your computer. It's for capturing video with software on your computer. Ezcap292 Medical Imaging Workstation supports UVC system, and work with most third video capture software well.

Record

Standalone mode - save video and photo in USB flash drive, USB HDD or 2.5 inch HDD directly

After connection done as above, screen should show the input video, and you also can get sound if you connect video output to TV or monitor. Please make sure the storage you select has been connected well, and LED next to the storage should turn on. ezcap292 Medical Imaging Workstation is ready now.

To record video is very simple --- press the VIDEO button (or REC on remote control), then you will see the red LED next to VIDEO button turns on. Press VIDEO button again (or STOP key on remote control) to stop recording, the LED turns flashes few seconds then turns off. Press PHOTO button (or PHOTO key on remote control), then you will see the red LED next to

VIDEO button turn on few seconds then turn off, it means device snapshot a photo and save in storage already. **It only works when NOT recording.**

Press UP/DOWN button to adjust the audio volume, and press LEFT/RIGHT button to adjust the microphone volume.

Press 720P/1080P/480P keys of remote control to select recording resolution and video output resolution.

Press WATERMARK key of remote control to add time watermark on recording file. Press again to cancel time watermark.


There are several reasons for no actions when press VIDEO button, you may need check:

- 1, The input video source disconnected.
- 2, The storage is out of memory.
- 3, The format for storage is not suitable, e.g. exFAT. ezcap292 support FAT, FAT32, and NTFS.
- 4, The storage is not connected well.
- 5, The storage is not ready, normally, it will take 4~8 seconds after plug it on.

Schedule recording:

Please set the SYSTEM TIME to the right local time before you set your schedule recording.

Press SCH REC of remote control, set the schedule time and others first, then click Select icon to display the task detail in task subject. Device will start and stop recording following your settings automatically.

Set the recording video quality & split:

Press MENU button (or MENU key on remote control), enter to following menu.


Clock: Set the time of system

Schedule Recording: Set the schedule record (same function as above)

Record Bitrate: Set recording video quality, set bit rate to High-16Mbps, Middle-10Mbps or Low-4Mbps.

Record File: Set record file to SPLIT or NO_SPLIT(NTFS). SPLIT means will cut the record file to 1.95GB each both in FAT and NTFS storage. NO_SPLIT(NTFS) means will save in one whole file in NTFS storage, but cut the record file to 1.95GB each in FAT storage.

Time WaterMark: Set time watermark with recorded file (same function as above)

please do not take the storage off when the LED is flashing. The recording will automatically stop if the USB flash out of space.

Sometimes you may find that some recorded video with problems, like: skip. It may because you don't click "safety remove" the device when disconnect from computer. You can try to format it to FAT32 if that problem keeps happen.

Video capture card mode – capture video with software on your computer

ezcap292 Medical Imaging Workstation also can work with a video capture card. It is USB 3.0 interface, and supports UVC system, can work with most of video capture and broadcast software well. You just connect USB 3.0 cable with PC port of back plate and free USB 3.0 port of computer, and run the application software, then choose capture device to "ezcap Medical HD" .

If video source includes audio, and get the bad audio quality, please set up the sample rate of recording device to 48KHz.

Connect device to computer, and right click **Speaker** icon in the bottom right corner of your computer screen, choose **Recording devices**. Double click **ezcap Medical HD**, at **Advanced**, choose the **sample rate** to **2 channel 16bit, 48000 Hz (DVD Quality)**.


Note: Most of all, you will get nice sound after you finish above settings. If you still get bad audio, it means the audio sample rate of your video source is 44.1KHz or 32KHz, please set here to 44.1 or 32 to make the sound without noise.

ezcap292 Medical Imaging Workstation does not only support Windows, but also supports Mac OS X and Linux, you just set the video capture device to "ezcap Medical HD" .


Playback / Check recorded file

Playback recorded file directly

After recording, you may want to playback the recorded file. Just press the PLAYBACK button (or PLAYBACK key on remote control), you will see the playback page on screen.


Choose PHOTO folder to check photo files, and choose MOVIE to check video files:

Use menu direction button (or direction buttons on remote control) to choose movie then press PLAY/PAUSE button to playback the recorded video.

You can press LAST, PLAY/PAUSE, NEXT, FAST BACKWARD, FAST FORWARD button to operate playback as usual.

Delete recorded video file: Use menu direction button (or keys of remote control) to choose the file you want to delete, press OK button to select it, then press DELETE button, press LEFT or RIGHT to choose Yes and press OK to delete the file.


Playback, check, or backup recorded file on your computer

You also can pull out USB flash drive or USB HDD from USB HOST, and connect to your computer to backup them, and playback with other video playback software.


And you can connect mini USB cable with READ HDD port and your computer, and **set device to USB mode (not HDD mode)**, then you can backup the video files in HDD to your computer, and playback with other video playback software.

other Settings


Press MENU button to enter menu operations, press direction keys to choose the settings. Set OSD Language.

Picture mode:

Sound mode:

Update firmware.

There are two firmware for this device – system part, record part.
Please make sure the device is on USB HOST mode before you update firmware.

1, Update firmware for SYSTEM PART: Copy the firmware file (MERGE.bin) to USB flash drive, plug USB drive to USB HOST, then press MENU button, use direction key to choose **Software Update (USB)** as above picture, then press OK and LEFT to update. Screen will display updating process from 1% to 100%, then screen turns black and display video again, it means device update firmware successfully.

2, Update firmware for RECORD PART: Copy the firmware file (jedi.img) in to USB flash drive, press POWER button to turn off device, plug USB drive to USB HOST, then press POWER button again to turn on device. The red indicator LED next to VIDEO button will flash few seconds, then LED turns off. It means the updating is finished. Plug USB flash drive to computer, and the jedi.img file disappears, it also means updating successfully.

Specifications

General

Recordable Devices	USB Flash Drive, USB Hard Drive, 2.5 inch Hard Drive. (supports FAT, FAT32 and NTFS format, and does not support exFAT format)
Display	4.3 inch LCD
Power supply	DC 12V3A
Operating temperature	41 to 104 deg F (5 to 40 deg C)
Storage temperature	-4 to 140 deg F (-20 to 60 deg C)

Operating humidity	5 至 80% (No condensation)
Dimention	262(L) X 172(W) X 120(H) mm
Weight	1.6 Kg
Accessories	Ezcap292 Medical Imaging Workstation X1, Remote Control X 1, IR Receiver cable X 1, USB 3.0 cable X 1, Mini USB cable X 1, Operation Manual X 1

Recording Format

Still Image	JPEG
Video	MPEG-4 H.264

Input/Output

Video Input	HDMI X 1, SDI(BNC) X 1, DVI X 1, VGA X 1, CVBS(BNC) X 1, YPbPr(BNC) X 1
Audio Input	HDMI(from HDMI video input) X 1, SDI(BNC)(from SDI video input) X 1, LINE IN(3.5mm stereo jack) X 1, MIC(3.5mm stereo jack) X 1
Video Output	HDMI X 1, *SDI(BNC) X 1, *CVBS(BNC) X 1 (Note: SDI OUT is only for SDI input, CVBS OUT is only for CVBS input)
Audio Output	HDMI(from HDMI video output) X 1, *SDI(BNC)(from SDI video output) X 1, LINE OUT(3.5mm stereo jack) X 1, (Note: SDI audio output is only for SDI input)
Other Interfaces	Video Foot Pedal port(3.5mm stereo jack) X 1, Photo Foot Pedal port(3.5mm stereo jack) X 1, IR Receiver cable port(3.5mm stereo jack) X 1, READ HDD port(Mini USB) X 1, PC port(USB

	3.0) X 1, Extension port(RS232) X 1, DC Power Supply port X 1
--	---

Recording Times for a 500GB Hard Drive (1920 X1080)

High quality (16Mbps)	75 hours
Medium quality (10Mbps)	115 hours
Low quality (5Mbps)	275 hours

*Features and specifications are subject to change without notice.

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