

THE HIGH-PERFORMANCE TORSIONAL ULTRASONIC SCALPEL



LOTUS – with torsional ultrasonic technology for enhanced safety

Ultrasonics - then and now

Ultrasonic scalpels have been used for laparoscopic surgery since the 1990s. In technological terms, an ultrasonic scalpel is a mechanical instrument with a vibrating blade.

Oscillating at a frequency of 36,000 Hz, it effectively acts as a simultaneous cutting and coagulation instrument.

Ultrasonic energy enables fast dissection without requiring the use of high frequency current.



The latest torsional ultrasonic technology

With the LOTUS system you benefit from the advantages of the latest generation of ultrasonic scalpels. The torsional energy generated in the LOTUS system reduces stray energy dissipation at the tip of the device, when compared to conventional longitudinal instruments

Vessels are sealed quickly and reliably with the LOTUS ultrasonic scalpel.

The patented torsional ultrasonic technology makes the LOTUS system especially efficient.

Enhanced Safety



Schematic depiction



Schematic depiction


Conventional ultrasonic instrument

Energy is fed in a linear direction through to the tip of the instrument. This leads to stray energy being dissapated at the tip of the instrument. Inadvertent distal penetration of tissue is possible.

LOTUS torsional ultrasonic scalpel

The LOTUS system's energy is perpendicular (90 degrees) to the axis of the blade. This coupled with the blade geometry focuses the energy into the jaw area. Lotus also offers enhanced safety by reducing the risk of inadvertent damage from stray distal energy.






LOTUS torsional transducer

Tool-less instrument assembly


Handpiece Simple Plug'n Play setup



The cutting speed of the LOTUS ultrasonic scalpel is directly related to the pressure on the trigger.

The high degree of tactile feedback enables precise cutting control by the surgeon.

Easy handling, precise operation



Efficient through and through

Dissecting shears

LOTUS dissecting shears are specifically designed for fast, precise haemostatic tissue dissection. The thin curved blade has focussed recesses and facilitates accurate dissection at the desired location.

Liver resector

The LOTUS liver resector is specifically made for use on liver parenchyma tissue. The larger contact surface creates a stronger haemostatic effect.



Dissecting shears

Narrow blade with focusing grooves for accurate dissection.



Liver resector

Larger contact surface for a stronger haemostatic effect.





LOTUS dissecting shears



^{*} only for LG4 Series 4 generators with software version Issue 6 or later. Previous versions must be used with the transducer CV3-400.

LOTUS liver resector





Optional accessories

Туре	REF
LOTUS LG4 foot switch	LF4
ARC CART, angeled equipment shelf for LOTUS (shelf for retrofitting to ARC CART)	902-924
LOTUS CART (device trolley), assembled Consisting of: 902-024 LOTUS CART, unassembled 902-921 Grip 902-912 Basket 902-100 Tool	902-070

Technical data

Technical data summary	LOTUS Series 4 generato

Line frequency	50/60 Hz
Mains voltage	100-240V (±10%)
Mains fuse	Internal mains fuse
Max. power consumption	150 VA
Width x height x depth	340 x 95 x 340mm
Weight	4.3 kg
Display	Thin-film technology
Classification according to EC Directive 93/42/EEC	IIb
Protection class according to EN 60601-1	I
Type according to EN 60601-1	BF
Ref	LG4
Output frequency	36 kHz
Output power	70 W (± 30 W)*
Duty type	Intermittent 3/30s (on/off)
	*

*depending on the transducer type

Ambient conditions for storage and transport

Temperature	-10°C to +50°C
Relative humidity	10 % to 90 %
Atmospheric pressure	500 hPa to 1060 hPa

Ambient conditions for operation

Temperature	+10°C to +30°C
Relative humidity	30 % to 75 %
Atmospheric pressure	810 hPa to 1060 hPa



BOWA-electronic GmbH & Co. KG Heinrich-Hertz-Strasse 4 – 10 72810 Gomaringen I Germany

Phone +49 (0) 7072-6002-0 Fax +49 (0) 7072-6002-33



