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EAR DROPS "DEWAX SQUALENE" INSTRUCTIONS FOR USE

Components.

Active components: hydrogen peroxide 35% – 0.86%, squalene – 0.1%.

Auxiliary components: glycerin, chlorhexidine digluconate.

Description.

A homogeneous liquid from light yellow to yellow in color, without foreign inclusions and impurities. Slight cloudiness or sediment is allowed.

Form of manufacture.

10 ml polymer dropper bottle, placed in a cardboard box along with the instructions for use. Not less than 280 doses.

Properties.

This product does not exert its effect through pharmacological, immunological, or metabolic processes, but rather through the mechanical and chemical action of its functional components.

The effect of "DEWAX SQUALENE" is due to the properties of its constituent components.

Hydrogen peroxide upon contact with organic substances (earwax), it breaks down into water and atomic (active) oxygen. This process is accompanied by the release of gas bubbles, creating a foaming effect. Active oxygen interacts with the components of the earwax plug, breaks down its structure, penetrates deep into the solid plug, breaks it down into smaller particles, and helps detach it from the walls of the ear canal.

Squalene acts as a cerumenolytic – a substance that softens and dissolves earwax, facilitating its removal by reducing the surface tension of the wax, allowing it to be expelled without additional mechanical manipulation. Squalene – a natural component of the hydrophilic layer that quickly forms a protective barrier on the skin's surface, preventing irritation and stimulating skin regeneration. A natural antioxidant, it possesses anti-inflammatory and antimicrobial properties.

Glycerin is used as a softening agent due to its properties as a substance that rapidly draws water from animal and plant tissues and as a substance that prevents drying out and freezing.

Chlorhexidine digluconate (preservative) is an antimicrobial component that significantly reduces the risk of microbial contamination that could cause product spoilage.

Ear drops "DEWAX SQUALENE" for external use.

Are used:

- for the prevention of earwax plugs formation;
- to dissolve earwax and maintain ear hygiene, including in cases of excessive earwax production;
- in case of itching or discomfort in the ear;
- for frequent users of headphones or earphones;
- to remove dust and moisture from the ear canal as a hygienic measure;
- for regular users of hearing aids;
- to protect the ears from chlorine compounds when swimming in a pool.

Method of application.

To remove earwax plugs in adults and children aged 6 and older, it is recommended to instill 2-4 drops into each ear canal 2-3 times a day for 3-4 days.

For daily hygiene and to prevent earwax buildup in adults and children aged 6 and older, it is recommended to instill 1 drop into each ear canal 2-3 times a month.

Precautions, warnings, and potential side effects.

To prevent discomfort resulting from contact between the ear canal skin and a cold solution, warm the bottle in your hands before use.

Usage by children is allowed only in the presence of adults.

In order to prevent spread of infection, one bottle should not be used by several people.

This product must not be used by people with individual hypersensitivity to its components.

Discontinue use at the first sign of a skin rash or other signs of local or systemic hypersensitivity.

Do not instill ear drops "DEWAX SQUALENE" into the eyes or nose.

Wait at least 15 minutes between using "DEWAX SQUALENE" and other medical devices or medicinal products.

Do not use this product if the eardrum is damaged or if damage is suspected. If the solution enters the middle ear, it may cause irritation, inflammation, and temporary hearing loss.

Excessive use may disrupt the natural protective balance of the external auditory canal, cause dryness, irritation, itching, or discomfort, and increase the risk of infectious complications.

Expiration date.

See on the packaging. After opening the bottle, it is recommended to use the product within 1 month.

Storage conditions.

Store at a temperature from +5 °C to +25 °C in a place protected from light and inaccessible to children.

Do not use after expiration of the date indicated at the package.

Manufacturer.

"Georg BioSystems" LLC.

3, Vokzalna str., Pobuzke village, Holovanivsk district, Kirovohrad region, 26555, Ukraine.

Address of manufacture.

64B, Dorozhna str., city of Kamianske, Dnipropetrovsk region, 51905, Ukraine.

Tel.: 0 800 501 161. E-mail: contact@georgbio.com

To prevent counterfeiting, the original packaging of the product contains an insert with a unique code in addition to the instructions for use.

Side 1	Side 2
	AA 00 00 00 00 00

THE MEANING OF THE SYMBOLS SHOWN ON THE MARKING:

	Use – by date
	Batch code
	Manufacturer
	Date of manufacture
	Temperature limit
	Consult instructions for use
	Mark of conformity with technical regulations
	Keep away from sun light