UA.TR.120

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EAR DROPS "OTOGUT SQUALENE" INSTRUCTIONS FOR USE

Components.

Active components: eucalyptus essential oil – 0.5%, calendula extract – 0.15%, sage extract – 0.1%, squalene – 0.1%.

Auxiliary components: glycerin, water, chlorhexidine digluconate, polyquaternium-37.

Description.

A homogeneous liquid from light yellow to yellow in color, without foreign inclusions and impurities. Slight cloudiness or sediment is allowed.

Form of manufacture.

10 ml polymer dropper bottle, placed in a cardboard box along with the instructions for use. Not less than 280 doses.

Properties.

Ear drops "OTOGUT SQUALENE" are a medical device intended for topical application in the external auditory canal and, as recommended by a doctor, in cases involving tissue care following procedures in the middle ear region.

The product's action is not achieved through pharmacological, immunological, or metabolic means, but rather through the physical and physicochemical properties of its components, which provide moisturization, softening, and coating of tissue surfaces, as well as protection against external irritants.

The product promotes the formation of a protective layer on the tissue surface, maintains the proper condition of the tissues of the external auditory canal, and creates conditions for natural tissue regeneration following irritation, water treatments, procedures, or surgical interventions in the ear area.

The product's components are selected based on their physicochemical properties, which ensure the care, moisturization, softening, and protection of ear tissues.

Eucalyptus essential oil inhibits the growth and reproduction of bacteria that can cause infectious otitis, reduces inflammation and swelling in the affected area, disinfects the application site, and prevents the spread of infection.

Calendula extract has antimicrobial properties. It exhibits antiviral activity against a wide range of microorganisms, including staphylococci and streptococci, which are common causes of ear infections. The flavonoids and organic acids in the extract help reduce swelling and hyperemia (redness) of the ear canal tissues, accelerate the regeneration of damaged skin and mucous membranes of the outer ear, and protect against the formation of rough keloids at the site of injury.

Sage extract has antiseptic properties, kills pathogenic bacteria-including staphylococci and streptococci, relieves irritation and swelling, and accelerates the healing of tissues and mucous membranes in the ear canal.

Squalene stimulates the regeneration of skin cells and mucous membranes. A natural antioxidant, it possesses antimicrobial properties. It has the ability to scavenge thermodynamically active singlet oxygen atoms formed as a result of oxidative reactions. In human skin cells, squalene absorbs free radicals and reduces damage caused by oxidative stress.

Glycerin is used as a softening agent due to its properties as a substance that rapidly extracts water from animal and plant tissues, and as a substance that prevents drying and freezing.

Chlorhexidine digluconate (preservative) is an antimicrobial component that significantly reduces the risk of microbial contamination that could cause product spoilage. The purpose of using preservatives is either to destroy bacteria, yeasts, and molds or to slow down their growth and development.

Indications for use.

Ear drops "OTOGUT SQUALENE" are used:

- for the care of the tissues of the external auditory canal in cases of irritation;
- to moisturize, soften, and protect the tissues of the external auditory canal after exposure to water, cold, and other external factors;
- to maintain the proper condition of the tissues of the external auditory canal following medical procedures;
- for caring for tissues following surgical procedures in the middle ear region – to create a protective environment and support natural tissue regeneration;
- as an adjunct to maintain the proper condition of ear tissues when experiencing dryness, irritation, or discomfort.

Method of application.

For children aged 6 years and older and adults instill 4 drops into the external auditory canal 2-3 times daily.

The duration of use should not exceed 10 days, unless otherwise recommended by a doctor.

When used following surgical procedures in the middle ear, the product should be used in accordance with the doctor's recommendations.

Precautions, warnings, and potential side effects.

To prevent discomfort caused by the cold solution coming into contact with the skin of the ear canal, warm the bottle in your hands before use. Use for children is permitted only under adult supervision.

Do not share a single bottle among multiple people.

Do not use this product if you have an individual hypersensitivity to any of its components.

If a skin rash, irritation, or other signs of local or systemic hypersensitivity occur, discontinue use.

Ear drops "OTOGUT SQUALENE" should not be instilled into the eyes or nose.

The interval between the use of this product and other medical devices or medications should be at least 15 minutes.

In cases of use following surgical procedures in the middle ear area, use the product as recommended by a doctor.

Expiration date.

See on the packaging. After opening the bottle, it is recommended to use the product within 5 weeks.

Storage conditions.

Store at a temperature from +5 °C to +25 °C in a place protected from light and inaccessible to children.

Do not use after expiration of the date indicated at the package.

Manufacturer.

"Georg BioSystems" LLC.

3, Vokzalna str., Pobuzke village, Holovanivsk district, Kirovohrad region, 26555, Ukraine.

Address of manufacture.

64B, Dorozhna str., city of Kamianske, Dnipropetrovsk region, 51905, Ukraine.

Tel.: 0 800 501 161. **E-mail:** contact@georgbio.com

To prevent counterfeiting, the original packaging of the product contains an insert with a unique code in addition to the instructions for use.

Side 1	Side 2
	AA 00 00 00 00 00

THE MEANING OF THE SYMBOLS SHOWN ON THE MARKING:

	Use – by date
	Batch code
	Manufacturer
	Date of manufacture
	Temperature limit
	Consult instructions for use
	Mark of conformity with technical regulations and identification number of the designated conformity assessment body
	Keep away from sun light