



SKINTACT ECG Electrode FS50

Aqua-Wet, Stress-foam backing material

KEY FEATURES

- foam backing materials protect sensor and gel from surgical and cleaning fluids
- the liquid gel is designed for fast pick up of the ECG signal even if inadequate skin prep is carried out
- all sensors have an high quality Ag/AgCl layer



PRODUCT

<i>Recomm. Application</i>	(short term) monitoring
<i>Characteristics</i>	disposable, pregelled, no latex, no PVC, non sterile
<i>Shelf life</i>	24 months unopened
<i>Storing conditions (min/max)</i>	+5° C/ +30° C
<i>X-ray translucent</i>	-

DIMENSIONS

<i>Electrode shape</i>	round – 50 mm
<i>Electrode size - max L/W [mm]</i>	-
<i>Total area [mm²]</i>	1.963
<i>Gel area [mm²]</i>	241
<i>Adhesive area [mm²]</i>	1.583

CLASSIFICATION AND STANDARDS

<i>Classification (MDD Art. 9)</i>	Class I
<i>Classification (CFR 21 870.2360.)</i>	Class II

MATERIALS

<i>Stud</i>	stainless steel
<i>Label</i>	PET-foil
<i>Backing material:</i>	PE-stress-foam
<i>Adhesive</i>	medical-grade acrylate
<i>Sensor</i>	silver/silverchloride (Ag/AgCl)
<i>Sponge</i>	PET-foam
<i>Gel</i>	Aqua-Wet (liquid gel)
<i>Release liner</i>	siliconized PET-foil (transparent)

MATERIALS PACKAGING

<i>Pouches, inner layer</i>	Polyethylene (PE)
<i>Pouches, centre layer</i>	Aluminium (Al)
<i>Pouches, outer layer</i>	Paper
<i>Boxes</i>	Cardboard

STANDARD PACKAGING

Item Number: 58028

<i>Pieces / card</i>	1
<i>Pieces / pouch</i>	30
<i>Pieces / box</i>	1200

Lot number and expiration date on every pouch and box

DIMENSIONS PACKAGING

<i>Box (L/W/H) [cm]</i>	39/ 29/ 22
-------------------------	------------

BIOCOMPATIBILITY

<i>Test - ISO 10993-1</i>	<i>Backing material (adhesive)</i>	<i>Gel</i>
<i>Cytotoxicity</i>	pass	pass
<i>Skin Irritation</i>	pass	pass
<i>Sensitization</i>	pass	pass

ELECTRICAL VALUES

	<i>Units</i>	<i>Typical Values</i>	<i>AAMI limits</i>
<i>DC-Offset</i>	[mV]	≤ 0,5	≤ 100
<i>DC Offset (5 sec after capacitor discharge)</i>	[mV]	≤ 10	≤ 100
<i>Recovery Slope</i>	[mV/sec]	≤ -0,2	± 1
<i>AC-Impedance with 10 Hertz</i>	[Ω]	≤ 100	≤ 2000
<i>Internal Noise Test</i>	[µV]	≤ 10	≤ 150
<i>Bias-Tolerance (over 24 h)</i>	[mV]	≤ 10	≤ 100